

可调无油真空泵

Adjustable Speed Oil-free Diaphragm Vacuum Pump



抽滤速度可调
The suction speed is adjustable



带灯开关，安全美观
Illuminated switch



过载、过热、短路保护
Overload,overheat and short circuit protection



真空表精度高，反应快
High precision vacuum gauge

技术参数

TECHNICAL DATA

型号 Model	VP-15S
类型 Type	负压 Negative pressure
抽气速率 Suction flow rate	0-15 L/min
极限真空度 Limit vacuum degree	0.08Mpa
噪音 Noise	< 45db
外壳保护等级 Enclosure ratings	IP21
电机功率 Power supply	30W
隔膜阀片 Diaphragm	进口橡胶 Rubber
电子保护 Protection devices	过载保护、过温保护、短路保护 Overload protection, over temperature protection, short circuit protection
毛重/净重 Gross/net weight	1.95/1.58 KG