


无油隔膜真空泵 - 30L

Oilless Diaphragm Vacuum Pump - 30L



正负压
Positive and negative pressure



带灯开关，安全美观
Illuminated switch



过载、过热、短路保护
Overload,overheat and short circuit protection



真空表精度高，反应快
High precision vacuum gauge

技术参数

TECHNICAL DATA

型号 Model	VP-30L
类型 Type	正负压 Positive and negative pressure
抽气速率 Suction flow rate	30 L/min
极限真空度 Limit vacuum degree	0.085Mpa
最大压力 Greatest pressure	5.0bar
噪音 Noise	< 50db
外壳保护等级 Enclosure ratings	IP21
电机功率 Power supply	50W
隔膜阀片 Diaphragm	进口橡胶 Rubber
电子保护 Protection devices	过载保护、过温保护、短路保护 Overload protection, over temperature protection, short circuit protection