

废液回收器

Liquid Vacuum Aspiration System



多种吸头适配

Work for muti purpose



数显调速旋钮

Digital display speed control knob



9档调速

9-speed speed regulation



可高温消毒

High temperature disinfection

技术参数

TECHNICAL DATA

| | | | |
|----------------------|--|--|------------|
| 型号 Model | VP-6Pro | 吸液速率 Fluid absorption rate | 1.6L / Min |
| 类型 Type | 负压抽气 Negative pressure suction | 极限真空度 Ultimate vacuum | 0.06Mpa |
| 控制方式 Control type | 数显旋钮·9档调速 Digital display knob·9 speeds | 吸液瓶直径 Diameter of collection bottle | 93mm |
| 抽气速率 Pumping rate | 6L / Min | 功率 Power | 10W |