

过滤装置

Vacuum Filtration Apparatus

技术参数

TECHNICAL DATA

接收瓶 Receiving flask	砂芯直径 Sand core diameter	砂芯孔径 Sand core aperture
250ml	20mm	10um
500ml	40mm	10um
1000ml	40mm	10um
2000ml	40mm	10um
5000ml	40mm	30-50um
10000ml	40mm	30-50um



高硼硅
Boro 3.3



密封性好
Tight sealing



适配GL-45试剂瓶
Suitable for bottles
with GL-45



大口径·抽滤迅速
Large diameter·Fast
suction



技术参数

TECHNICAL DATA

接收瓶 Receiving flask	砂芯直径 Sand core diameter	滤膜 Filter membrane	砂芯孔径 Sand core aperture
2000ml	80mm	90-100mm	10um
5000ml	80mm	90-100mm	10um
10000ml	80mm	90-100mm	10um