

水浴锅

Water Bath



数显恒温
Digital display thermostatic



控温精度高
High precise temperature control



不锈钢一次成型内胆
Stainless steel tank



抗蚀性强
High corrosion resistance

技术参数

TECHNICAL DATA

型号 Model	WB100-1 WB100-1F	WB100-2 WB100-2F	WB100-4 WB100-4F	WB100-6 WB100-6F
温度范围 Temperature range	室温 - 99.9°C Room temp - 99.9°C			
功率 Power	300W	600W	800W	1200W
精度 Accuracy	0.1°C			
工作室尺寸 Working chamber size	16*16*13cm	30*15*14cm	30*30*14cm	45*30*13cm

* WB100-F系列为防干烧系列 / WB100-F series are anti-dry series