


## ELECTRIC STIRRER

- ◆Cantilever type design, can prevent the liquid into the machine
- ◆Large range of cooling the body surface, uniform rapid cooling, prevent the operation temperature rises
- ◆Motivation is big, is applicable to work long hours
- ◆Especially suitable for the mixing of large capacity with high viscosity fluid used in laboratory.



• GZ-20

• GZ-20A



• HS-160C



- ◆Microcomputer intelligent control, LED display speed,time
- ◆Avoid maintenance, low noise, DC brushless motor drive
- ◆Transmission stirrer, it is convenient to change the container during the experiment
- ◆Automatic overload protection, guarantee the safe operation of the experimental high load for a long time
- ◆Microcomputer control, timing control state, the end of the work after five minutes, can automatically cut off the machine power supply, safe, energy-saving.
- ◆Strong stirring, suitable for high viscosity medium and paste work for a long hours

Model	GZ-20	GZ-20A
Speed adjustable range	60-500rpm 240-2000rpm	60-500rpm 240-2000rpm
Power	120W	120W
Holding the agitator	0.5-10mm	0.5-10mm
Display of speed	----	YES
Voltage	110V 60Hz / 220V 50Hz	110V 60Hz / 220V 50Hz
Dimension(mm)	300x220x700	300x220x700
Package(mm)	430x380x200	430x380x200
G.W./N.W.	12/10kg	12/10kg

Model	HS-160C	HS-200C
Speed range	50-2000rpm	50-2000rpm
Power	160W	200W
Max torque	100Nm	150Nm
Timer	0-999min	0-999min
Max stirring viscosity	20000Ps	30000Ps
Maxmixing capacity	20L	30L
Clamping stirring rod	∅ 0.5-10mm	∅ 0.5-10mm
Voltage	110V 60Hz / 220V 50Hz	110V 60Hz / 220V 50Hz
Dimension(mm)	300x220x700	300x220x700
Package(mm)	430x380x140	430x380x140
G.W./N.W.	12/10kg	14/12kg



• JJ-1



• JJ-5



• ES-300C

Product name	Electric Stirrer	Speed measuring Electric Stirrer	LCD constant speed Electric Stirrer
Model	JJ-1	JJ-5	ES-300C
Power	60W 100W 160W 200W	60W 100W 160W 200W	300W
Speed range	Start up-3000rpm	Start up-3000rpm	Start up-2000rpm
Timer	0-120min or all the time	0-120min or all the time	0-999min or all the time
Voltage	110V 60Hz/220V 50Hz	110V 60Hz/220V 50Hz	110V 60Hz/220V 50Hz
Max stirrer water capacity	20L	20L	50L
Characteristic	----	Digital display speed	Drill chuck, ∅ 0.5-10mm Brushless motor drive
Dimension(mm)	370x200x640	370x200x640	370x340x680
Package(mm)	460x330x170	480x430x170	460x380x170
G.W./N.W.	8/7kg	10/9kg	13/12kg



• A-88



• FSH-2A



• JJ-3

Product name	Tissue Triturator	High Speed Homogenizer	Six-in-one Electric Stirrer
Model	A-88 (JJ-2)	FSH-2A	JJ-3
Power	200W	185W	6x60W
Capacity	1000ml	5ml,10ml,15ml,20ml	6x5000ml
Speed range	start up ~12000 rpm	Blade No.1: 8000~18000 rpm Blade No.2: 8000~20000 rpm	start up ~2800 rpm
Dimension(mm)	230x230x420	300x250x400	1000x260x650
Package(mm)	280x280x470	350x300x430	1100x350x700
G.W./N.W.	8/7kg	8/7kg	45/34kg