

# WATER BATH

- ◆Microprocessor control , touch panel control
- ◆Multi purpose: room temperature to 99°C
- ◆Optional function: room temperature to 200 °C
- ◆One-timer-molding cistern , easy to clean
- ◆LCD display temperature setpoint and temperature measurements, with 0.1 °C resolution
- ◆Inside and outside pump cycle,inside circulation ensures temperature uniformity
- ◆With power-off protection,low water level protection and alarm
- ◆With the higher and lower temperature alarm setting function
- ◆With temperature run away protection and alarm function



• ZEN-TW



• 601



- ◆Microprocessor control,
- ◆Digital display setting /inner temperature
- ◆Good insulation performance with stainless steel cover
- ◆With Inner and outer cycle
- ◆Suitable for thermostatic , baking, dewaxing, sterilization

Model	ZEN-TW	601
Capacity	27L	20L
Inner size(mm)	460x300x200	400x280x180
Working size(mm)	340x300x200	300x280x180
Power	1200W	1000W
Temperature range	RT+5--99°C	RT+5 -99°C
Temperature accuracy	± 0.5°C	± 0.5°C
Timer	0-999min	0-999min
Pump flow	8-10L/min	6-8L/min
Voltage	110V 60Hz / 220V 50Hz	110V 60Hz /220V 50Hz
Dimension(mm)	530x330x430	480x380x400
Package(mm)	600x390x490	530x440x460
G.W./N.W.	16/14kg	14/12kg

- ◆Microprocessor PID control,
- ◆Digital display setting and actual temperature
- ◆Digital display running timer
- ◆PID self-tuning , over temperature alarm
- ◆Automatic alarm after operation
- ◆Stainless steel circle cover, easy to operation
- ◆Inside and outside pump cycle



• DHC-1010



• DLC-1120



- ◆Microcomputer control,touch panel control
- ◆LCD display temperature setpoint and temperature measurements, with 0.1 °C resolution
- ◆Inside and outside pump cycle
- ◆With power-off ,low water level protection and alarm
- ◆Over current protection, high temperature and low temperature alarm function, empty operated protection
- ◆With temperature runaway protection and alarm, parameter memory ,temperature correction function

Model	DHC-1005	DHC-1010	DHC-1020
Capacity	17L	17L	17L
Working size(mm)	300x260x220	300x260x220	300x260x220
Diameter (mm)	φ215	φ215	φ215
Heating power	1000W	1000W	1000W
Power	200W	260W	400W
Temperature range	-5--99°C	-10--99°C	-20--99°C
Temperature homogeneity	± 0.5°C	± 0.5°C	± 0.5°C
Voltage	220V 50Hz	220V 50Hz	220V 50Hz
Dimension(mm)	450x390x860	450x390x860	450x390x860
Package(mm)	580x520x910	580x520x910	580x520x910
G.W./N.W.	54/42kg	58/46kg	62/50kg

Model	DLC-1120	DLC-0620
Capacity	20L	14L
Inner size(mm)	400x320x160	350x270x150
Open size(mm)	160x160	160x160
power	1500W	1000W
Temperature range	-20--99°C	-20--99°C
Temperature accuracy	± 0.1°C	± 0.1°C
Pump circulation mode	Inner and outer cycle	Inner and outer cycle
Pump flow	8-10L/min	8-10L/min
Voltage	110V 60Hz / 220V 50Hz	110V 60Hz / 220V 50Hz
Dimension(mm)	500x430x900	450x380x820
Package(mm)	620x500x950	550x460x870
G.W./N.W.	58/45kg	49/37kg