

WATER BATH

- ◆Microprocessor control,
- ◆Digital display inner temperature
- ◆Stainless steel interior
- ◆The aluminum ring cover
- ◆Suitable for thermostatic , baking, dewaxing, sterilization



• HH-S4

- ◆Microprocessor control,
- ◆Digital display setting and actual temperature
- ◆Digital display running timer
- ◆PID self-tuning , over temperature alarm
- ◆Automatic alarm after operation
- ◆Suitable for thermostatic , baking, dewaxing, sterilization



• DK-S26



• WH-4



- ◆Microprocessor control,
- ◆Digital display inner temperature
- ◆One-timer-molding cistern,
- ◆The aluminum ring cover
- ◆Stainless steel roof cover (optional)
- ◆Suitable for thermostatic , baking, dewaxing, sterilization



• DK-8A

• DK-600B

- ◆Microprocessor control,
- ◆Digital display setting and actual temperature and time
- ◆PID self-tuning , over temperature alarm
- ◆Automatic alarm after operation
- ◆Good insulation performance with stainless steel roof cover
- ◆With bakelite handle
- ◆DK-600B with double layer cover.

Model	HH-S1	HH-S2	HH-S4	HH-S6	HH-S8
Capacity	3L	6L	10L	16L	22L
Power	300W	500W	1000W	1500W	2000W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1.5℃	±1.5℃
Temperature range	RT+5 -99℃				
Temperature homogeneity	±1℃				
Voltage	110V 60Hz /220V 50Hz				
Working size(mm)	160x160x120	300x160x120	300x300x120	460x300x120	600x300x120
Dimension(mm)	230x230x240	370x230x240	370x370x240	530x370x240	670x370x240
Package(mm)	290x290x280	410x280x280	420x410x280	610x430x280	730x430x280
G.W./N.W.	3.5/3kg	4.5/4kg	6/5 kg	7.5/6.5kg	9/8kg

Model	WH-1	WH-2	WH-4	WH-6
Capacity	3L	6L	12L	16L
Power	200W	500W	1000W	1500W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1.5℃
Temperature range	RT+5 -99℃			
Temperature homogeneity	±1℃			
Voltage	110V 60Hz /220V 50Hz			
Working size(mm)	150x140x120	300x150x120	320x320x120	460x300x120
Dimension(mm)	180x180x240	330x180x240	360x330x240	540x330x250
Package(mm)	230x220x300	380x230x280	410x400x280	590x400x290
G.W./N.W.	3/2.5kg	4/3kg	6/5 kg	8/7 kg

Model	DK-S12	DK-S24	DK-S26	DK-S28
Capacity	6L	10L	16L	22L
Power	500W	1000W	1500W	2000W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1℃
Temperature range	RT+5 -99℃			
Temperature homogeneity	±1℃			
Timer	0-999min			
Voltage	110V 60Hz /220V 50Hz			
Working size(mm)	300x160x120	300x300x120	460x300x120	600x300x120
Dimension(mm)	470x230x240	470x370x240	630x370x240	770x370x240
Package(mm)	540x290x280	540x430x280	710x430x280	880x430x280
G.W./N.W.	5.5/5kg	7/6.5 kg	11/10 kg	13/12 kg

Model	DK-320	DK-420	DK-8A	DK-600B
Capacity	6L	9L	23L	20L
Power	500W	1000W	1200W	1500W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1℃
Temperature range	RT+5 -99℃			
Temperature homogeneity	±0.5℃			
Timer	0-999min			
Voltage	110V 60Hz /220V 50Hz			
Working size(mm)	300x160x120	400x180x130	450x300x170	580x280x130
Dimension(mm)	470x230x260	570x270x280	630x370x320	760x360x260
Package(mm)	530x300x310	630x320x320	670x430x360	810x420x320
G.W./N.W.	7/6kg	8/7kg	12/11kg	15/14kg