

**Product Characteristics**


**hot plate**



DB-2



DB-3A



HP-500A

Product name	Hot plate			
Model	DB-1	DB-2	DB-3	DB-4
Heating power	500W	1000W	2000W	2500W
Platform Dimension(mm)	200x150mm	300x200mm	380x200mm	400x300mm
Voltage	220V 50Hz	220V 50Hz	220V 50Hz	220V 50Hz
Dimension(mm)	200x150x140/4kg	300x200x140/6kg	380x200x140/8kg	400x300x140/12.5kg

Product name	Intelligent temperature hot plate			
Model	DB-1A	DB-2A	DB-3A	DB-4A
Temperature range	RT~300℃	RT~300℃	RT~300℃	RT~300℃
Heating power	500W	1000W	2000W	2500W
Platform Dimension(mm)	200x150mm	300x200mm	380x200mm	400x300mm
Voltage	220V 50Hz	220V 50Hz	220V 50Hz	220V 50Hz
Dimension(mm)	300x150x140/5kg	400x200x140/7kg	480x200x140/9.5kg	500x300x140/14kg

Product name	Ceramic hot plate		
Model	HP-500A	HP-1000A	HP-1500A
Max Temperature	250℃	250℃	250℃
Heating power	500W	1000W	1500W
Platform Dimension(mm)	184x184mm	300x200mm	300x300mm
Voltage	220V 50Hz	220V 50Hz	220V 50Hz
Dimension(mm)	300x190x180/ 6kg	300x300x180/8kg	450x300x180/14kg