

CENTRIFUGE

- ◆ Strong Plastic case with metal lid
- ◆ With mechanical door lock
- ◆ With safety switch, Stop working when open the door
- ◆ Scope of supply includes the rotor
- ◆ High speed for short centrifugation timers
- ◆ Small footprint
- ◆ LED display speed, Stepless speed adjustment



• LC-04A

- ◆ Strong Plastic case with metal lid
- ◆ High Brushless motor, with mechanical door lock
- ◆ Over speed and imbalance protection
- ◆ With safety switch, stop working when open the door
- ◆ Stepless speed adjustment, LED display speed and timer
- ◆ Error code display, speed up and down adjustment



• LC-04C



• LC-04B

- ◆ Strong Plastic case with metal lid
- ◆ With safety switch, Stop working when open the door
- ◆ Pointer show speed
- ◆ Stepless speed adjustment
- ◆ Mechanical timer
- ◆ Scope of supply includes the rotor
- ◆ High speed for short centrifugation timers



• LC-04L

- ◆ Strong Plastic case
- ◆ Brushless motor, LCD display speed, RCF and timer
- ◆ Over speed alarm and imbalance protection
- ◆ Touch panel control, automatic stop once lid open
- ◆ Error code display, Accel /Decel rates adjustment
- ◆ Electronic Door Lock, Equipped with emergency locking device
- ◆ Suit for qualitative analysis of blood serum and plasma

Model	LC-04A	LC-04B
Max RPM	4000 RPM, LED display speed	4000 RPM, Pointer show speed
Max RCF	2325xg	2325xg
Timer	0-30min	0-30min
Rotor capacity	20mlx12 (standard)	20mlx12 (standard)
Voltage	110V 60Hz /220V50Hz	110V 60Hz /220V50Hz
Dimension	280x310x265mm	280x310x265mm
Package	340x350x330mm	340x350x330mm
G.W./N.W.	10/9kg	10/9kg

Model	LC-04C	LC-04L
Max RPM	4000 RPM	4000 RPM
Max RCF	2325xg	2325xg
Timer	0-99min	0-99min
Rotor capacity	20mlx12 (standard)	20mlx12 (standard)
Voltage	110V 60Hz /220V 50Hz	110V 60Hz /220V 50Hz
Dimension	280x310x265mm	320x340x280mm
Package	340x350x330mm	410x430x350mm
G.W./N.W.	10/9kg	11/10kg

Rotor NO.	Rotor capacity	Max RCF	LC-04A	LC-04B
A 042012 (standard)	20mlx12	2325xg	◆	◆
A 042006	20mlx6	1790xg	*	*
A042008	20mlx8	2325xg	*	*
A040518	15mlx18	2325xg	*	*
A041012	10mlx12	2325xg	*	*

Rotor NO.	Rotor capacity	Max RCF	LC-04C	LC-04L
A 042012 (standard)	20mlx12	2325xg	◆	◆
A 042006	20mlx6	1790xg	*	*
A042008	20mlx8	2325xg	*	*
A040518	15mlx18	2325xg	*	*
A041012	10mlx12	2325xg	*	*
A040512	5ml eppendorf tube x12	1790xg	*	*
A040812	8ml x12	2325xg	*	*