

CENTRIFUGE

- ◆ High strength plastic shelf , simple structure and simple operation.
- ◆ With mechanical door lock and safety switch, stop working when open the door.
- ◆ Pointer show speed, stepless speed adjustment
- ◆ Viewing area can view the running statement.



• LC-04B Plus

- ◆ High strength plastic lid, small volume, light weight
- ◆ brushless motor drive
- ◆ Large screen backlight LCD shows speed/running time/Accel / Decel
- ◆ With RCF/RPM toggle key, review the RCF at any time
- ◆ Error code , speed Accel/ Decel rates 10 programs
- ◆ Advanced wide voltage design, RPM is not affected by voltage fluctuation, high precision
- ◆ Over speed alarm, finished alarm
- ◆ Mechanical door lock and safety switch, stop open the door



• LC-04C Plus



Model	LC-04B Plus	LC-04C Plus
Max RPM	4000 RPM	4000 RPM
Max RCF	2325xg	2325xg
Timer	0-30min	0-99min
Rotor capacity	20mlx12 (Standard)	20mlx12 (Standard)
Voltage	110V 60Hz /220V 50Hz	100-240V 50-60Hz
Dimension	320x340x280mm	320x340x280mm
Package	410x430x350mm	410x430x350mm
G.W./N.W.	9/8kg	8/7kg

Rotor NO.	Rotor Capacity	Max RCF	LC-04B Plus	LC-04C Plus
A 042012	20mlx12	2325xg	◆	◆
A042008	20mlx8	2325xg	*	*
A040518	5mlx18	2325xg	*	*
A041012	10mlx12	2325xg	*	*
A 042006	20mlx6	1790xg	*	*
A.040812	8mlx12	2325xg	---	*

PRP Centrifuge

- ◆ Strong Plastic case with microprocessor control
- ◆ LED display speed and timer
- ◆ DC Brushless motor,
- ◆ Over speed and finished alarm
- ◆ Mechanical door lock ,the cover open and will stop
- ◆ Error code display, speed up and down adjustment
- ◆ Used for PRP/PRF,
- ◆ Suit for centrifuge the 10 ml syringe directly
- ◆ It is dedicated product for the beauty salon



• LC-04P



• LC-04P plus



- ◆ Strong plastic case with microprocessor control
- ◆ LCD display speed/RCF / time /Acceleration / Decel/ rotor NO. /Door station
- ◆ DC brushless motor, Over speed alarm and finished alarm
- ◆ mechanical door lock ,the cover open and will stop
- ◆ 10 groups of programs with memory function.
- ◆ Logo display when turn on the power, Programs can be programmed according to user parameters
- ◆ Suit for centrifuge the 10ml syringe directly
- ◆ For the beauty salon

Model	LC-04P	LC-04P plus
Max RPM	4000 RPM	4000 RPM
Max RCF	1780xg	1780xg
Timer	0-99min	0-99min59s
Rotor capacity	Angle rotor 15ml /10mlx8	Angle rotor 15ml /10mlx8
Voltage	100-240V 50-60Hz	100-240V 50-60Hz
Dimension	260x300x205mm	260x300x205mm
Package	330x350x280mm	330x350x280mm
G.W./N.W.	6/5kg	6/5kg