

# CENTRIFUGE



• LC-04W

- ◆Microcomputer control, brushless motor variable frequency
- ◆Large screen backlight LCD shows
- ◆10 groups of ACC /DEC rate to choose and memory
- ◆Electronic door lock control,with emergency unlock device
- ◆Stainless steel chamber
- ◆Imbalanced protection, self locking, overspeed protection, speed no signal protection, motor holzer fault protection, high or low power supply voltage protection, over-current protection

Model	LC-04W
Max RPM	4000 RPM
Max RCF	3150xg
Timer	0-99min
Max Rotor capacity	400ml
Voltage	110V 60Hz /220V 50Hz
Dimension	410x480x340mm
Package	510x580x460mm
G.W./N.W.	39/31kg

Description	Rotor NO.	Rotor Capacity	Max RPM/RCF
<b>LC-04W</b>	S041524	15mlx24	4000rpm 2860xg
Optical Swing Rotor 1	S041516	15mlx16	4000rpm 2860xg
Optical Swing Rotor 2	S041508	15mlx8	4000rpm 2610xg
Optical Swing Rotor 3	S041506	15mlx6	4000rpm 2540xg
Optical Swing Rotor 4	S0410004	100mlx4	4000rpm 3150xg
Optical Swing Rotor 5	S045008	50mlx8	4000rpm 2860xg
Optical Swing Rotor 6	S045004	50mlx4	4000rpm 2860xg
Optical Swing Rotor 7	S043006	30mlx6	4000rpm 2860xg
Optical Swing Rotor 8	S041024	10mlx24	4000rpm 2500xg
Optical Swing Rotor 9	S040804	8mlx6	4000rpm 2500xg
Optical Swing Rotor 10	S040712	7mlx12	4000rpm 2500xg

## PRP Centrifuge

- ◆Microprocessor control ,LCD display speed and time ,RCF.
- ◆Brushless frequency motor, free maintenance, quick in speed up and down.
- ◆Used for fat transplantation purification
- ◆Suit for screw syringe of 50ml, 20ml and 10ml.



• LC-04M

- ◆Strong metal case, big capacity, low noise
- ◆Swing rotor, suit for 4 different tubes
- ◆Stainless steel chamber and hood
- ◆High power brushless motor,
- ◆LED display speed and timer, stepless speed adjustment
- ◆Over speed and imbalance protection
- ◆With safety switch and electromagnetic door locks
- ◆Error code display, speed up and down adjustment



• TDL-4C/TDL-5A

Model	LC-04M
Max RPM	5000RPM
Max RCF	4400xg
Timer	0-99min
Swing capacity	50ml Syringex4 (standard) /3500rpm 15ml x16 (optical) /4200rpm
Voltage	110V 60Hz / 220V 50Hz
Dimension	450x590x350mm
Package	550x670x460mm
G.W./N.W.	48/36kg

Model	TDL-4C	TDL-5A
Max RPM	4000 RPM	5000 RPM
Max RCF	2755xg	4390xg
Timer	0-99min	0-99min
Rotor capacity	Swing 100mlx4 ,50mlx8, 15mlx16,10mlx24	Swing 250mlx4 ,50mlx8, 15mlx28,10mlx36
Voltage	110V 60Hz /220V50Hz	110V 60Hz /220V50Hz
Dimension	450x550x450mm	450x550x450mm
Package	540x660x560mm	540x660x560mm
G.W./N.W.	58/44kg	62/48kg