

CENTRIFUGE

- ◆ Strong Plastic case with metal lid
- ◆ With mechanical door lock and safety switch
- ◆ Scope of supply includes the rotor
- ◆ High speed for short centrifugation timers
- ◆ Small footprint
- ◆ Stepless speed adjustment
- ◆ LED display speed
- ◆ Suit for qualitative analysis of blood serum and plasma



• LC-05A

- ◆ Strong Plastic case with metal lid
- ◆ High power brushless motor
- ◆ Over speed and imbalance alarm
- ◆ Stepless speed adjustment , LED display speed and timer
- ◆ With mechanical door lock and safety switch
- ◆ Error code display, speed up and down adjustment
- ◆ Suit for qualitative analysis of blood serum and plasma



• LC-05C



• LC-05B



- ◆ Strong Plastic case with metal lid
- ◆ With mechanical door lock and safety switch
- ◆ Scope of supply includes the rotor
- ◆ High speed for short centrifugation timers
- ◆ Stepless speed adjustment
- ◆ Suit for qualitative analysis of blood serum and plasma



• LC-06C

- ◆ Strong Plastic case with metal lid
- ◆ High power brushless motor
- ◆ Over speed and imbalance alarm
- ◆ Stepless speed adjustment , LED display speed and timer
- ◆ With mechanical door lock and safety switch
- ◆ Error code display, speed up and down adjustment
- ◆ Suit for qualitative analysis of blood serum and plasma

Model	LC-05A	LC-05B
Max RPM	5000 RPM	5000 RPM
Max RCF	3600xg	3600xg
Timer	0-30min	0-30min
Rotor capacity	20mlx12 (standard)	20mlx12 (standard)
Voltage	110V 60Hz /220V 50Hz	110V 60Hz /220V 50Hz
Dimension	280x310x265mm	280x310x265mm
Package	340x350x330mm	340x350x330mm
G.W./N.W.	10/9kg	10/9kg

Rotor NO.	Rotor capacity	Max RCF	LC-05A	LC-05B
A 052012 (standard)	20mlx12	3600xg	◆	◆
A 052006	20mlx6	2760xg	*	*
A052008	20mlx8	3600xg	*	*
A051518	15mlx18	3600xg	*	*
A051012	10mlx12	3600xg	*	*

Model	LC-05C	LC-06C
Max RPM	5000 RPM	6000 RPM
Max RCF	3600xg	5200xg
Timer	0-99min	0-99min
Rotor capacity	20mlx8 (standard)	15mlx12 (standard)
Voltage	110V 60Hz /220V 50Hz	110V 60Hz /220V 50Hz
Dimension	280x310x265mm	280x310x265mm
Package	340x350x330mm	340x350x330mm
G.W./N.W.	10/9kg	11/10kg

Description	Rotor No.	Rotor capacity	Max RCF
LC-05C 5000rpm	A052008 (standard)	20mlx8	3600xg
	A 052006	20mlx6	2760xg
	A 052012	20mlx12	3600xg
	A051518	15mlx18	3600xg
	A051012	10mlx12	3600xg
LC-06C 6000rpm	A 061512 (standard)	15mlx12	5200xg
	A 062006	20mlx6	4000xg
	A062008	20mlx8	5200xg
	A061518	15mlx18	5200xg
	A051012	10mlx12	5200xg