

AUTOMATIC STEEL STERILIZER

- ◆ Hand wheel translation type quick opening structure
- ◆ High quality stainless steel material
- ◆ Microcomputer automatic control, arbitrary setting sterilization parameters
- ◆ Safety interlocking device
- ◆ Digital display and Touch type button
- ◆ Standardized test interface
- ◆ Self expanding seal
- ◆ With automatic protection function: over temperature protection;
- ◆ Over pressure self discharge protection; low water level protection;
- ◆ Anti dry burning
- ◆ With stainless steel basket
- ◆ Sterilization end buzzer reminder, automatic stop
- ◆ Automatic discharge of cold air, automatic exhaust steam after sterilization
- ◆ Optional external printer



• LDZX-50B

Model	LDZX-40B LDZX-40BD LDZX-40BP	LDZX-50B LDZX-50BD LDZX-50BP	LDZX-75B LDZX-75BD LDZX-75BP	LDZX-100B LDZX-100BD LDZX-100BP	LDZX-150B LDZX-150BD LDZX-150BP
Capacity	40L	50L	75L	100L	150L
Sterilization room size	ø350×445mm	ø350×525mm	ø400×625mm	ø450×650mm	ø500×760mm
Maximum pressure	0.22Mpa				
Maximum temperature	134℃				
Maximum timer	9999min				
Heating power	3.5KW	3.5KW	4KW	4.5KW	6.0KW
Power supply	220V 50Hz	220V 50Hz	220V 50Hz	220V 50Hz	380V 50Hz
Basket	ø320×400×1	ø320×240×2	ø360×280×2	ø410×280×2	ø450×240×3
Dimension(mm)	550×540×1000	550×540×1100	580×570×1100	630×620×1100	630×620×1460
Package(mm)	630×580×1130	630×580×1230	650×610×1350	690×660×1400	690×660×1760
G.W./N.W.	79/58kg	88/64kg	95/72kg	110/83kg	140/105kg

* Remark: LDZX-BD serial : with drying function.
LDZX-BP serial : FO Value printing function with sterilization effect.

Automatic Table Top Steam Sterilizer

- ◆ Quality stainless steel liner, sealed cover
- ◆ Microcomputer control, automatic cycle, arbitrary sterilization procedures, Touch buttons
- ◆ Safety interlocking device, Open top water tank
- ◆ Quipped with standard test interface
- ◆ Europe CLASS N.standard
- ◆ Over temperature , over voltage , low water level protection, anti dry burning
- ◆ Sterilization end buzzer reminder, automatic stop
- ◆ Automatic cooling air, sterilization end automatic exhaust steam
- ◆ Built-in steam generator to provide saturated steam quickly



• B18C



• B45C-V

- ◆ Automatic door switch
- ◆ Equipped with water quality test
- ◆ Microcomputer control, Safety interlocking device
- ◆ Equipped with standard test interface
- ◆ Open type top water tank
- ◆ European class B standard
- ◆ With automatic protection: over temperature , over pressure, low water level protection, anti dry burning , safety door lock.
- ◆ Three times per-vacuum
- ◆ Sterilization end buzzer reminder, automatic stop
- ◆ Automatic discharge of cold air, automatic exhaust steam after sterilization.
- ◆ Built-in steam generator to provide saturated steam quickly
- * Optional external printer
- * The optional USB function and U disk

Model	B8C	B12C	B16C	B18C	B23C
Capacity	8L	12L	16L	18L	23L
Power	1500W	1800W	1800W	2200W	2200W
Power supply	110V 60Hz /220V 50Hz				
Max pressure	0.22 MPa				
Working temperature	134 °C				
Sterilization room size(mm)	ø170×360	ø200×385	ø230×390	ø250×370	ø250×470
Dimension(mm)	410×530×400	410×530×400	480×550×450	480×560×450	480×560×450
Package(mm)	660×500×450	660×500×450	680×570×510	680×570×510	770×560×500
G.W./N.W.	43/38kg	45/40kg	52/44kg	56/48kg	61/53kg

Model	B8C-V	B16C-V	B18C-V	B23C-V	B45C-V
Capacity	8L	16L	18L	23L	45L
Power	1500W	1800W	2200W	2200W	3500W
Power supply	110V 60Hz /220V 50Hz				
Max pressure	0.22 MPa				
Working temperature	134 °C				
Sterilization room size(mm)	ø170×360	ø230×390	ø250×370	ø250×470	ø315×615
Dimension(mm)	410×530×400	480×550×450	480×550×450	480×560×450	700×790×750
Package(mm)	660×500×450	680×570×510	680×570×510	770×560×500	1000×780×790
G.W./N.W.	47/42kg	55/47kg	60/52kg	65/58kg	134/105kg