

DRY OVEN

Dry Oven

- ◆With hidden door locks , silent fans
- ◆Mechanical overheat protection
- ◆Electricity leakage protection.
- ◆With double glass window
- ◆Microprocessor control,
- ◆Digital display setting and actual temperature and time
- ◆Over temperature alarm , finished alarm



• OV-9053A

Dry /Incubator (Dual-purpose)

- ◆Microprocessor control, oven/incubator hierarchical control
- ◆With double glass window
- ◆Digital display setting and actual temperature
- ◆Digital display running timer
- ◆PID self-tuning, over temperature alarm
- ◆Automatic alarm after operation
- ◆Stainless-steel interior, chamber corners are semicircle



• PPH-030A

Model	OV-9023A	OV-9053A	OV-9073A	OV-9123A
Capacity	25L	50L	70L	105L
Working size(mm)	300x300x275	415x375x345	450x350x450	550x370x520
Power	500W	1000W	1200W	1500W
Temperature range	RT+10-200°C			
Temperature resolution	0.1°C			
Temperature fluctuation	±1°C			
Timer	0-9999min			
Voltage	110V 60Hz /220V 50Hz			
Dimension(mm)	430x500x620	510x650x690	480x660x800	690x630x860
Package(mm)	530x600x720	610x750x820	580x760x930	730x790x990
G.W./N.W.	42/30kg	59/49kg	67/55kg	84/65kg

Model	PPH-030A	PPH-050A	PPH-070A	PPH-140A	PPH-240A
Capacity	40L	50L	70L	140L	225L
Working size(mm)	350x350x350	415x375x345	450x350x450	550x450x550	500x600x750
Power	850W	850W	1100W	1600W	2100W
Temperature range	As incubator: RT+5~80°C, As oven 80~250°C				
Temperature resolution	0.1°C				
Temperature fluctuation	As incubator: ±0.5°C, As oven: ±1°C				
Timer	0-9999min				
Voltage	110V 60Hz /220V 50Hz				
Dimension(mm)	620x630x520	690x610x530	620x710x620	720x780x750	770x840x980
Package(mm)	700x710x620	790x660x620	720x780x730	830x880x860	940x870x1140
G.W./N.W.	54/38kg	58/45kg	64/56kg	88/72kg	116/90kg

Dry Oven

- ◆Microprocessor control,
- ◆With double glass window
- ◆Digital display setting and actual temperature
- ◆Digital display running timer
- ◆PID self-tuning, over temperature alarm
- ◆Automatic alarm after operation
- ◆Stainless-steel interior, chamber corners are semicircle



• DHG-9053A

Model	DHG-9023A	DHG-9053A	DHG-9123A	DHG-9423A
Capacity	25L	50L	105L	420L
Working size(mm)	300x300x275	415x375x345	550x370x520	600x550x1300
Power	500W	1000W	1500W	5000W
Temperature range	50~200°C			
Temperature resolution	0.1°C			
Temperature fluctuation	±1°C			
Timer	0-999min			
Voltage	110V 60Hz /220V 50Hz			
Dimension(mm)	580x540x460	690x600x540	825x660x690	740x830x1650
Package(mm)	670x610x530	770x670x660	920x730x820	850x920x1750
G.W./N.W.	41/29kg	56/40kg	78/58kg	210/170kg

Model	DHG-9030A	DHG-9070A	DHG-9140A	DHG-9240A	DHG-9624A
Capacity	40L	70L	140L	225L	640L
Working size(mm)	350x350x350	450x350x450	550x450x550	500x600x750	800x800x1000
Power	500W	1000W	1500W	2000W	5000W
Temperature range	50~300°C				
Temperature resolution	0.1°C				
Temperature fluctuation	±1°C				
Timer	0-999min				
Voltage	110V 60Hz /220V 50Hz				
Dimension(mm)	620x630x520	620x680x620	720x780x750	770x840x980	940x1090x1380
Package(mm)	680x720x610	720x750x740	830x860x870	900x870x1130	1000x1190x1510
G.W./N.W.	52/38kg	62/48kg	86/64kg	111/89kg	227/200kg