



• DZF-6020/6050



• DZF-6090/6210

- ◆Microprocessor control,
- ◆With glass window door,
- ◆Digital display setting and actual temperature
- ◆Digital display running timer
- ◆PID self-tuning, over temperature alarm
- ◆Automatic alarm after operation
- ◆Suitable for drying, baking, dewaxing, sterilization, heat treatment under vacuum conditions
- ◆DZF-6090/6210 with vacuum pump

Model	DZF-6020	DZF-6050	DZF-6090	DZF-6210
Capacity	25L	50L	90L	210L
Working size(mm)	300x300x275	415x375x345	450x450x450	560x640x600
Power	1000W	1400W	1600W	2160W
Vacuum pump	Without	Without	Yes	Yes
Temperature range	50~180℃	50~200℃	50~250℃	50~250℃
Temperature fluctuation	± 1%			
Vacuum degree	<133Pa			
Timer	0-999min			
Voltage	220V 50Hz			
Dimension(mm)	590x460x470	710x540x550	610x650x1430	810x760x1590
Package(mm)	670x540x570	810x620x650	750x740x1580	890x840x1680
G.W./N.W.	54/42kg	76/60kg	163/117kg	195/140kg

Shaking Incubator

- ◆Microprocessor control temperature and timer
- ◆LED display temperature ,speed and timer
- ◆Stainless steel inner, with door lock.
- ◆stepless shaking speed adjustable
- ◆The glass door is convenient for observation.
- ◆The interior lighting is convenient for observation
- ◆Over temperature alarm, Finished work alarm
- ◆BOD culture, seed-breeding test
- ◆Thermostatic microorganism, plant cultivation



• BS-1E



• SHTR-140

- ◆ Microprocessor control temperature and timer
- ◆ LCD display temperature ,speed and timer
- ◆ Stainless steel inner,
- ◆ Stepless shaking speed adjustable ,brushless motor driver
- ◆ The glass window is convenient for observation.
- ◆ Over temperature alarm, Finished work alarm
- ◆ With power module failure, motor stall, Hall error, bus voltage undervoltage, overvoltage, communication failure alarm function

Model	BS-1E	BS-2F
Capacity	140L	280L
Temperature range	5-50℃	5-50℃
Temperature accuracy	± 1℃	± 1℃
Shaking speed	Start up ---300rpm	Start up---300rpm
Amplitude	20mm(Orbital)	20mm(Orbital)
Bottles capacity	Tube ø16x200,100mlx15,200mlx9	(Tube ø16x200,100mlx15,200mlx9)x2
Heating power	300W	400W
Compressor power	130W	180W
Timer	0-9999min	0-9999min
Voltage	220V 50Hz	220V 50Hz
Dimension(mm)	570x590x1100	1140x600x1100
Package(mm)	640x700x1310	1230x700x1200
G.W./N.W.	90/62kg	160/110kg

Model	SHTR-140	SHTR-280
Capacity	140L	280L
Temperature range	4-60℃	4-60℃
Temperature accuracy	± 0.5℃	± 0.5℃
Shaking speed	Start up ---300rpm	Start up---300rpm
Amplitude	20mm(Orbital)	20mm(Orbital)
Bottles capacity	Tube ø16x200,100mlx15,200mlx9	(Tube ø16x200,100mlx15,200mlx9)x2
Heating power	300W	400W
Compressor power	200W	400W
Timer	0-9999min	0-9999min
Voltage	220V 50Hz	220V 50Hz
Dimension(mm)	600x630x1150	1180x630x1150
Package(mm)	700x730x1250	1280x730x1250
G.W./N.W.	105/75kg	190/135kg