

WATER BATH

- ◆Microprocessor control,
- ◆Digital display inner temperature
- ◆Good insulation performance with stainless steel roof cover
- ◆With bakelite handle



• HH-W420



• WB-18S

- ◆Microcomputer PID intelligent control and LCD backlight display temperature and time
- ◆The upper cover is equipped with ring butterfly cover
- ◆Stainless steel inner gallbladder tensile forming process, the inner wall is smooth, the product is more beautiful
- ◆Over temperature protection, leakage current protection, low water alarm protection
- ◆Prevent dry burning heating components, effectively protect heating elements
- ◆Optional stainless steel flat cover or roof cover

Model	HH-W420	HH-W600
Capacity	9L	20L
Power	1000W	1500W
Temperature range	RT+5 -99℃	RT+5 -99℃
Temperature fluctuation	±1℃	±1℃
Temperature homogeneity	±1℃	±1℃
Voltage	110V 60Hz /220V 50Hz	110V 60Hz /220V 50Hz
Working size(mm)	400x180x130	580x280x130
Dimension(mm)	480x260x280	660x360x280
Package(mm)	530x320x320	710x420x330
G.W./N.W.	7/6kg	11/10kg

Model	WB-5S	WB-12S	WB-18S
Volume	5L	12L	18L
Working Size (mm)	300x150x120	320x320x120	460x300x130
Power	500W	1000W	1500W
Temperature range	RT+5-99℃	RT+5-99℃	RT+5-99℃
Temperature accuracy	±0.3℃	±0.3℃	±0.3℃
Timer	0-9999Hour	0-9999Hour	0-9999Hour
Voltage	220V 50Hz	220V 50Hz	220V 50Hz
Dimension(mm)	350x200x270	370x370x270	530x320x290
Package(mm)	400x260x320	430x430x320	590x380x350
G.W./N.W.	6/5kg	10/9kg	13/12kg

Cycle Water Bath

- ◆Microprocessor control
- ◆Digital display inner temperature
- ◆One-timer-molding cistern
- ◆Good insulation performance with stainless steel cover
- ◆With inner cycle
- ◆Prevent dry burning heating components effectively protect heating elements
- ◆Optional stainless steel flat cover or roof cover



• WH-6C



• HH-4

Oil Bath

- ◆Microprocessor PID control
- ◆Digital display setting and actual temperature
- ◆Digital display running timer
- ◆PID self-tuning , over temperature alarm
- ◆Automatic alarm after operation
- ◆Annular butterfly cover, easy to operate
- ◆High efficiency heat preservation material, small heat loss



Model	WH-4C	WH-6C	WH-60C
Capacity	12L	16L	62L
Power	1000W	1500W	2500W
Pump power	40W	40W	60W
Circulation pump flow	≥4L/min	≥4L/min	≥8-10L/min
Temperature fluctuation	±0.1℃	±0.1℃	±0.5℃
Temperature range	RT+5 -99℃		
Temperature homogeneity	±0.5℃		
Voltage	110V 60Hz /220V 50Hz		
Working size(mm)	320x320x120	460x300x120	620x500x200
Dimension(mm)	460x330x240	640x330x250	700x450x320
Package(mm)	530x390x280	710x390x290	800x550x370
G.W./N.W.	9/8kg	12/11kg	49/35kg

Model	HH-1	HH-2	HH-4	HH-6	HH-8
Capacity	3L	6L	10L	16L	22L
Working size(mm)	160x160x120	300x160x120	300x300x120	460x300x120	600x300x120
Power	500W	1000W	1500W	2000W	2500W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1.5℃	±1.5℃
Temperature range	50-300℃				
Temperature homogeneity	±1℃				
Voltage	110V 60Hz /220V 50Hz				
Dimension(mm)	330x230x240	470x230x240	470x370x240	630x370x240	770x370x240
Package(mm)	410x280x290	530x290x280	530x430x280	690x430x280	840x450x280
G.W./N.W.	6/5.5kg	9/8kg	11/10kg	14/13kg	17/16kg