

INCUBATOR

Anaerobic Incubator

- ◆Microprocessor control
- ◆Digital display setting and actual temperature
- ◆Combine of operation room, sample room, Gas circuit, circuit control system and Soluble wax sterilizer.
- ◆With UV-lamps , with special gloves
- ◆Gas circuit can adjust flow and input various gas needed.
- ◆Stainless steel interior
- ◆The front window made of transparency and Impact resistance of special glass
- ◆Operation room is equipped with special vaccination rods sterilizer, soluble wax disinfection device eliminate oxygen catalys



• YQX-II



• GNP-9050



Water Jacket Incubator

- ◆Water jacket, double stainless steel interior
- ◆Over temperature alarm
- ◆Water level high or low, automatic alarm
- ◆Microprocessor control
- ◆Digital display setting and actual temperature
- ◆Digital display running timer, PID self-tuning
- ◆Automatic alarm after operation

Laboratory Incubator

- ◆Microprocessor control,
- ◆Digital display setting and actual temperature
- ◆Over temperature alarm, PID self-tuning
- ◆Digital display running timer,
- ◆Stainless steel interior
- ◆Electricity leakage protection.
- ◆With glass inner door, With hidden door locks,
- ◆Mechanical overheat protection, cut off the heating sources
- ◆Automatic alarm after operation



• IB-9052A



• DNP-9052A



- ◆Microprocessor control
- ◆Digital display setting and actual temperature
- ◆Over temperature alarm, PID self-tuning
- ◆Digital display running timer,
- ◆Stainless steel interior
- ◆With glass observe window
- ◆Automatic alarm after operation

Model	YQX-II
Working size (mm)	800x600x650
Chamber size (mm)	250x190x290
Power	950W
Temperature range	RT+5~50°C
Temperature resolution	±1°C
Time to Anaerobic state of sample room	<5min
Time to Anaerobic state of operation room	<1hour
Anaerobic state hold time	More than 12hours if stop providing trace gas mixture for operation room
Voltage	220V 50Hz
Dimension(mm)	1200x750x1330
Package(mm)	1300x850x1480
G.W./N.W.	230/174kg

Model	GNP-9050	GNP-9080	GNP-9160	GNP-9270
Capacity	50L	80L	160L	270L
Working size (mm)	350x350x410	400x400x500	500x500x650	600x600x750
Power	320W	470W	620W	770W
Temperature range	RT+5~65°C			
Temperature resolution	0.1°C			
Temperature fluctuation	±0.5°C			
Timer	0-999min			
Voltage	110V 60Hz /220V 50Hz			
Dimension(mm)	520x480x750	570x530x840	680x640x980	780x730x1080
Package(mm)	650x550x830	670x620x950	780x730x1120	880x820x1200
G.W./N.W.	53/42kg	58/45kg	78/59kg	100/79kg

Model	IB-9025A	IB-9052A	IB-9082A	IB-9162A	IB-9272A
Capacity	15L	50L	80L	160L	270L
Working size (mm)	250x250x250	350x350x410	400x400x500	500x500x650	600x600x750
Power	150W	250W	300W	400W	600W
Temperature range	RT+5~65°C				
Temperature resolution	0.1°C				
Temperature fluctuation	±0.5°C				
Timer	0-999min				
Voltage	110V 60Hz /220V 50Hz				
Dimension(mm)	380x370x580	480x470x740	530x530x840	640x630x1000	730x720x1090
Package(mm)	490x470x640	580x580x800	630x650x980	730x710x1120	880x860x1280
G.W./N.W.	35/21kg	46/35kg	62/47kg	90/64kg	112/86kg

Model	DNP-9025A	DNP-9030A	DNP-9052A	DNP-9082A	DNP-9162A	DNP-9272A
Capacity	15L	30L	50L	80L	160L	270L
Working size (mm)	250x250x250	300x300x350	350x350x410	400x400x500	500x500x650	600x600x750
Power	150W	200W	250W	300W	400W	600W
Temperature range	RT+5~65°C					
Temperature resolution	0.1°C					
Temperature fluctuation	±0.5°C					
Timer	0-999min					
Voltage	110V 60Hz /220V 50Hz					
Dimension(mm)	370x380x540	420x430x640	470x490x700	520x540x800	620x640x940	720x740x1010
Package(mm)	460x470x610	510x520x720	570x580x770	620x630x900	720x730x1080	810x820x1170
G.W./N.W.	30/20kg	38/28kg	48/36kg	59/44kg	77/60kg	103/82kg