

A

NO.	ROTOR NO.	PICTURE	CAPACITY
1	A042006		20mlx6
2	A042012		20mlx12
3	A041508		15mlx8
4	A041518		15mlx18
5	A041012		10mlx12
6	A042024		20mlx24
7	A045006		50mlx6
8	A160212		1.5/2mlx12
9	A160506		5mlx6
10	A160012		0.5mlx12
11	A130506		1.5mlx12 and 5mlx6
12	A130012		1.5mlx12 and 0.5mlx12
13	A120224		1.5/2mlx24
14	A130508		5mlx8
15	A121012		10mlx12
16	A130048		0.5mlx48
17	A131536		1.5mlx36
18	A1010004		100mlx4
19	A115006		50mlx6

S





NO.	ROTOR NO.	PICTURE	CAPACITY
20	S0410004		100mlx4 15mlx16 50mlx8 10mlx24
21	S0525004		250mlx4 15mlx28 50mlx8 10mlx36
22	S0510004		100mlx4
23	S055004		50mlx4
24	S052008		20mlx8
25	S042006		20mlx6
26	S121524		24 capillary tube
27	S042012		20mlx12
28	S041516		15mlx16
29	S040796		7mlx96
30	S041524		15mlx24
31	S041532		15mlx32
32	S045008		50mlx8
33	S040548		5mlx48
34	S0410008		100mlx8
35	S0450006		500mlx6
36	S0480004		800mlx4
37	S0405148		5mlx148
38	S044802		48 Microplatex2

§ And sample of our products in the continuous improvement of product structure, appearance, parameters and configuration may differ, please understand.



Test tubes Dimension (mm)

Centrifuge Adapter

Picture							
Capacity	800ml	500ml	250ml	100ml	50ml	20ml	15ml
Inner diameter mm	100	90	69.5	42.3	30.8	19/18.5	13
Length mm	106	100	100	96	91	101/93	92
Max insert size(mm)	∅ 100×100	∅ 90×105	∅ 69×110	∅ 42×105	∅ 30×103	∅ 18.5×113 ∅ 18×110	∅ 12.6×100

Micro-Centrifuge Tube

Picture							
Capacity	10ml	7ml	5ml	5ml	2ml	1.5ml	0.5ml
Outer diameter mm	15.5	12.9	13.3	16.5	10.3	10.3	7.5
Length mm	81	59	53	60	41	41	32

Screw Cap Centrifuge tube

Vacuum Tube

Picture								
Capacity	100ml	50ml	15ml	10ml	8ml	7ml	5ml	4ml
Outer diameter mm	32.5	28.6	16.5	16	15.5	16	13	13
Length mm	113	113	118	102	130	100	100	75