

Electrophoresis Power Supply



HT-600C



HT-600D

Universal Power Supply

Input range: 110-250V, suitable for use in different countries

Four sets of output jacks: High current can connect four electrophoresis tanks simultaneously

Type of Output: constant voltage, constant current, and constant power, continuously adjustable

Output range: HT-600C: 5-600V 1-600mA 5-300W

HT-600D: 5-600V 1-1200mA 5-500W

Increment: Voltage (1V) Current (1mA) Power (1W)

Timing range: 1 minute to 99 hours and 59 minutes

Temperature detection: -20 °C~+85 °C

Step-by-step control function: editable 10 groups, each group of up to make 10 program automatic connection operation Protection function: when overcurrent, overload, over temperature, no load, short circuit

Micro current: After timed operation, unmanned shutdown can maintain micro current output

Display screen: 7-inch large touch screen, automatic generation of Numeric keypad during setting, display timing, temperature, operation status and various prompt operations

Automatic setting: When a certain indicator of constant voltage, current, and power is set, the other two items automatically increase to the maximum value, No need for manual adjustment.

External dimensions: 30 * 30 * 13cm

