

Sequencing electrophoresis-Cell



HT-CX01

Suitable for SSR, AFLP, SSCP, HA and other experiments
Pure flat aluminum plate treated by special process, dissipate heat evenly, band without "smile" phenomenon
Large sample volume, 96 PCR samples can be electrophoresed at one time, improving work efficiency
High-quality red shark tooth comb, easy to observe and load samples, not easy to change samples
The buffer outlet of the upper tank and the movable bottom tank, easy to remove buffer
The upper and lower tank covers are designed to open the cover and cut off the power, ensure the safety of the experiment
99.99% high-purity platinum electrode, the best conductivity

Technical Specification:
Glass size: 35*35cm
Sample comb: 68、100wells
Buffer volume: 750mL

Name	Unit	Qty
Body tank	pc	1
Gel glass concave	pc	2
Gel glass plate	pc	2
Shark tooth comb 68-well	pc	2
Sharktooth comb 100-well	pc	2
Spacer	pc	4
Cable	pc	1



HT-CX02

Using buffer cooling method, the gel temperature is uniform, and the strips are neat, no deformation
Large sample volume, 96 PCR samples can be electrophoresed at one time, improving work efficiency
High-quality red shark tooth comb, easy to observe and load the sample, no sample change
There are buffer outlets in the upper and lower tanks, which is convenient for cleaning the buffer
Open the cover and cut off the power to ensure the safety of the experiment
99.99% high-purity platinum electrode, the best conductivity and corrosion resistance

Technical Specification:
Glass size (W*L) : 38*32cm
Sample comb: 68、100wells
Buffer volume: 1500mL

Standard configuration

Name	Unit	Qty
Body tank	pc	1
Gel glass concave	pc	2
Gel glass plate	pc	2
Shark tooth comb 68-well	pc	2
Shark tooth comb 100-well	pc	2
Spacer	pc	4
Cable	pc	1