

# Electrophoresis Power Supply



**HT-ECP3000**

## High voltage electrophoresis power supply

- This product is suitable for large and medium-sized horizontal and vertical electrophoresis, sequencing, bidirectional
- Commonly used electrophoresis such as isoelectric focusing.
- Output type: constant voltage, constant current, constant power, continuously adjustable
- Display: LCA display screen with backlight (128 \* 64 pixels)
- Microprocessor intelligent control has:
  - Overload, no-load, leakage and other protective functions
  - Can store 10 commonly used electrophoresis programs
  - Automatic memory function
  - Standard and timed operation function
  - Automatic judgment function
  - Constant voltage, constant current, and constant power intelligent prompts

## Technical parameter

Output range: 20-3000V, 1-200mA,  
1-200W

Resolution: voltage 1V, current 1mA, power 1W

Timing range: 1 minute to 99 hours and 59 minutes