



Vertical Electrophoresis



HT-Scz02

vertical electrophoresis tank (double plates)



HT-Scz03

vertical electrophoresis tank (four plates)

Applicable to gel electrophoresis of various proteins and SSR molecular markers
It can electrophoresis 2-4pieces medium-sized gel at the same time, meet the needs of high-throughput electrophoresis

Thickened gel glass plate, solid and durable

Integrated design of concave plate and edge strip

Special sealing mold to ensure no glue leakageThe buffer discharge outlet of the upper and lower tanks, easy to remove residual liquid

No additional electric field interference on both sides, with straight stripes

99.99% high-purity platinum electrode with optimal conductivity

Technical Parameters:

Glass plate size: 18×18cm

Gel size : 15×15cm

Gel thickness: 1.0mm

Sample comb: 20、34、50 (optional, add samples for multi-channel pipettes)

Standard configuration (HT-Scz02)

Name	Unit	Qty
Body tank	pc	1
Gel glass concave	pc	2
Gel glass plate	pc	2
20-well comb	pc	2
34-well comb	pc	2
Drainage pipe	pc	2
cables	pc	1
Sealing frame	pc	2
Stationery folder	pc	8

HT-Scz03 standard configuration, except for the main body of the electrophoresis tank, drainage pipes, and cables, twice the standard configuration of HT-Scz02