

Es5000 Intelligent Electronic Balance

After blood collection, blood bags can not be completely the same amount, ES5000 intelligent balancing comes to support the use of centrifuge, accurate and efficient configuration, to ensure the quality of blood separation and to improve the life of the centrifuge.



Main Features

- Large capacity of the tray easy for set the sample
- Easy to operate and reading, microcomputer control (ES5000-6S upgrade to PC control), high measurement accuracy.
- One machine with two functions, not only can measure also can be weighed
- Dual screen display makes readings easier and more convenient. (only for ES5000-6S models)
- Selection of silicone rubber, non-toxic and tasteless including: 2×3g、2×4g、1×5g、1×10g、1×20g.

Specification

Model	ES5000-2	ES5000-4	ES5000-6	ES5000-6S
Single weight fo tray	5000g	5000g	5000g	5000g
Number of tray	2	4	6	6
Readable resolution	1g	1g	1g	1g
Repeatable	1g	1g	1g	1g
Linear	1g	1g	1g	1g
Ambient temp	10-35℃	10-35℃	10-35℃	10-35℃
Sensitivity	0.05	0.05	0.05	0.05
Stable time	1.5second	1.5second	1.5second	1.5second
Single tray dimension	175×175mm	175×175mm	175×175mm	175×175mm
Power supply	AC220V 50Hz	AC220V 50Hz	AC220V 50Hz	AC220V50Hz
Dimendion	365×440×110mm	550×400×110mm	735×400×110mm	735×400×110mm