

Food Safety Centrifuge

Specification

Model	H/T85MM
Max Speed	8500r/min
Max Capacity	6×50mL
Max RCF	7500×g
Speed Adjustment	±50rrpm steps(can be customized)
Timer Range	1min-99h59min, with continuous run function
Control and drive system	Large torque AC variable frequency motor, microcomputer control
Total Power	400W
Power supply	AC220V 50Hz
Weight	35kg
Dimention (LxWxH)	710×630×350mm
Package dimensions	



Rotors Specification

Figure	Parameter
 6×50mL	Max Speed : 8500r/min Max RCF : 7500×g Serial No : H001-003-008-301 Bio-isolated sealed rotor cover