

Beauty Centrifuge

Used in plastic surgery, dermatology and other fields. It can be used for PRP injection transplantation. It uses 10mL, 20mL, 50mL conventional syringes, specially designed for cosmetologists and plastic surgeons in hospitals ,beauty salons ,health centers.

Specification

Model	TD5A Beauty Centrifuge
Max Speed	4000r/min
Max RCF	2100×g
Max Capacity	4×50mL
Speed Adjustment	±50rrpm steps (can be customized)
Timer Range	1min-99h59min, with continuous run function
Control and drive system	Large torque AC variable frequency motor, microcomputer control
Total Power	500W
Power supply	AC220V 50Hz
Noise level	≤60dB (A)
Weight	35kg
Dimention (LxWxH)	600×540×360mm
Package dimensions	

Special rotors and adapters can be customized according to user requirements.

