

# Cell smear centrifuge

Cell Smear Centrifuge is used to take erythrocyte-serum test, identify antigens and antibodies, and give other relative tests in immune hematology examination and research laboratories. also it can be used for cell purification by washing. This centrifuge machine is widely applicable for urology, cytology, microbiology, hematology, immunohistochemistry, clinical examination and so on.

## Specification

Model	TDZ4 Cell smear centrifuge
Max Speed	4000r/min
Capacity	12Cards
Max RCF	1824×g
Speed Adjustment	±50rrpm steps (can be customized)
Timer Range	1min-99h59min , with continuous run function
Control and drive system	Large torque AC variable frequency motor, microcomputer control
Power supply	AC220V 50Hz
Total Power	120w
Weight	20Kg
Dimendion (LxWxH)	500×450×280mm
Package dimensions	



## Rotors Specification

Figure



细胞涂片转子 12位

转子参数

Max Speed : 4000r/min  
Max RCF : 1824×g  
Serial No : H010-003-077-301