

Low Speed Urine Sediment Centrifuge

Low-speed multi-pipe rack urine sediment centrifuge is mainly used as urine sediment analysis instrument ,widely used in clinical medicine ,biochemistry,immunology and other fields

Specification

Model	TDZ5 Urine Sediment Centrifuge
Max Speed	4000r/min
Max Capacity	32×15mL
Max RCF	2960×g
Speed Adjustment	± 50rrpm steps (can be customized)
Timer Range	1min-99h59min, with continuous run function
Control and drive system	Large torque AC variable frequency motor, microcomputer control
Total Power	450W
Power supply	AC220V 50Hz
Weight	35kg
Dimention (LxWxH)	600×540×360mm

可根据用户要求，定制各种专用转子及adapter。

Rotors Specification

Figure	Parameter
 32 × 15mL	Max Speed : 4000r/min Max RCF : 2960×g Serial No : H008-002-065-201 Cup serial No : H002-002-050-208