

# Bench-top Dairy Centrifuge



Dairy centrifuge is the special equipment for milk and dairy standard for analytical methods. Applicable to pasteurization, cover Bolaw to separate fat milk and dairy products, measurement and analysis. Having the function of heating thermostat, and the temperature was maintained at 50 °C or more, which not only can reduce the time of the sample water bath, as well as the test results more reliable. Mainly applicable in dairy enterprises for testing, is commonly used in dairy inspection equipment. With TFT true-color LCD wide-

## Specification

Model	TLW5R	TLW5
Rotor type	Angle rotor	
Number of tubes	8 tubes	
Max Speed	1500r/min	
Max RCF	420×g	
Speed Adjustment	±50rpm steps (can be customized)	
Timer Range	1min-99h59min, with continuous run function	
Control and drive system	Large torque AC variable frequency motor, microcomputer control	
Temp control range	Room temp ~ +65°C	None
Total Power	1.5kW	400W
Noise level	≤60dB (A)	
Power supply	AC220V 50Hz	
Weight	35kg	
Dimension (L×W×H)	520×415×330mm	
Serial No	H019-003-100-301	
Package dimensions		