

K30



K30 series Dry Bath Incubator an economic type heating thermostat, equipped with a circular heating block.

Technical Data

Model	K30
Temperature Control Range	RT+5°C~100°C
Timing Range	1min~99h59min
Temperature Control Accuracy	±0.5°C
Heating Time(20°C to 100°C)	≤15min
Certificate	CE
Dimension(L x W x H) mm	225 x 160 x 197

Net Weight	3kg
------------	-----

Accessories

Code	Type	Description
AS-01031-01	A-30	96×0.2ml tubes
AS-01031-02	B-30	29×0.5ml+32×1.5ml tubes
AS-01031-03	C-30	61×0.5ml tubes
AS-01031-04	D-30	41×1.5ml tubes
AS-01031-05	E-30	41×2.0ml tubes
AS-01031-06	F-30	61xΦ6mm tubes
AS-01031-07	H-30	28xΦ12mm tubes
AS-01031-08	I-30	28xΦ13mm tubes
AS-01031-09	J-30	21xΦ15mm tubes
AS-01031-10	K-30	21xΦ16mm tubes
AS-01031-11	L-30	17xΦ18mm tubes
AS-01031-12	M-30	14xΦ20mm tubes
AS-01031-13	N-30	9xΦ25mm tubes
AS-01031-14	O-30	4xΦ32mm tubes