

Table Top Steam Sterilizer



Use

The table type sterilizers are safely auto-controlled, which are designed for the dental clinics, hospitals, school lab, scientific research institutions and laboratory. And mostly suitable for sterilizing of surgical, dental and eye instruments, glassware medicine, culturemedium and biological dressing, food, etc.

Characteristics

- ◆ Sterilizing course: Automatic sterilization controlled by computer, easy to operate.
 - ◆ 134°C, suitable for 4–6 minutes flash sterilizing.
 - ◆ Safe protection of water lacking.
- ◆ With over heat and over pressure protection.
 - ◆ After sterilization it will automatically cut the power and then alarm.
 - ◆ Fully stainless steel structure.

Table Top Steam Sterilizer



Technical Data

MODEL	TM-T16J	TM-T24J
Technical Data		
Sterilizing chamber volume	16L φ 250 × 360mm	24L φ 270 × 410mm
Design pressure	0.3MPa	
Design temperature	150°C	
Max. working pressure	0.25MPa	
Max. working temperature	138°C	
Timer range	0~60min	
Heat average	≤ ± 1°C	
Power	2KW/AC220V.50Hz	
Sterilizing plate(mm)	340 × 200 × 30(3piece)	
Dimension(mm)	545 × 328 × 423	560 × 350 × 500
Exterior package dimension(mm)	620 × 410 × 570	620 × 410 × 570
G.W/N.W	17kg / 15kg	20kg / 18kg