

## Table Top Steam Sterilizer



### Use

The table type sterilizers are safely auto-controlled, which are designed for the dental clinics, hospitals, school lab, scientific research institutions and laboratory. And mostly suitable for sterilizing of surgical, dental and eye instruments, glassware medicine, culture medium and biological dressing, food, etc.

### Characteristics

- ◆ Sterilizing course: Automatic sterilization controlled by computer, easy to operate.
- ◆ Maximum temperature: 134°C, suitable for 4-6 minutes flash sterilizing.
- ◆ Safe protection of water lacking.
- ◆ With over heat and over pressure protection.
- ◆ After sterilization it will automatically cut the power and then alarm.
- ◆ Fully stainless steel structure.

## Table Top Steam Sterilizer



### Technical Data

MODEL	TM-T16J	TM-T24J
Technical Data		
Sterilizing chamber volume	16L φ250×360mm	24L φ270×410mm
Working pressure	0.22MPa	
Working temperature	134°C	
Temperature adjust range	105~134°C	
Timer range	0~60min	
Heat average	≤ ± 1°C	
Power	1.5KW / AC220V.50Hz	
Sterilizing plate(mm)	340×200×30(3piece)	
Dimension(mm)	545×328×423	560×350×500
Exterior package dimension(mm)	620×410×570	620×410×570
G.W/N.W	17kg / 15kg	20kg / 18kg