

Portable Pressure Steam Sterilizer



Use

The product utilizes saturated steam to rapidly and efficiently sterilize articles. It can be used in hospitals, public health centers, health stations, clinics in factories and mines, scientific research institutes, etc, to sterilize medical equipment, surgical dressings, glass utensils, solutions, substrata. It can also be used as cooker in plateau areas and be used to make high quality drinking water in industrial and mining enterprises.

Characteristics

- ◆ Fully stainless steel structure.
 - ◆ Hand wheel type of quick-open door structure.
 - ◆ LED display of working status, touch type key.
- ◆ Automatically shut off with beep reminding after sterilization.
 - ◆ Electric heated.
 - ◆ Easy to operate, safe and reliable.

LED display automation



Portable Pressure Steam Sterilizer

Technical Data

MODEL		
Technical Data	YX-280E18	YX-280E24
Sterilizing volume	18L φ 280 × 260mm	24L φ 280 × 390mm
Working pressure	0.14~0.16MPa	
Working temperature	126℃	
Temperature adjust range	105~126℃	
Timer range	0~99min	
Max. Safety pressure	0.165MPa	
Heat average	≤ ± 1℃	
Power	AC220V.50Hz / 2KW	
Transport dimension(mm)	420 × 420 × 630	420 × 450 × 750
G.W/N.W	14kg / 12kg	15kg / 14kg