

Bante MS400C

Hotplate Magnetic Stirrer



Features

- Ceramic coated top plate is easy to clean.
- Large LCD display clearly shows the running time, temperature and stirring speed.
- Push-button speed control prevents the stir bar from accidentally stalling.
- External temperature sensor controls the heater to either heat up or maintain the temperature at the set value.
- Single point temperature calibration.

Stirrer Includes

- TP-100P temperature probe
- Stir bar
- Stand set
- Fuse

Specifications

Model	MS400C
Maximum Stirring Volume	0~2000 mL (H ₂ O)
Stirring Speed	0~2500 rpm
Timer Range	0 to 999 minutes
Heater	300W
Hotplate Temperature	Max. 300°C
Top Plate Size	135(L)×135(W)mm
Top Plate Material	Ceramic coating
Display	Custom LCD (95×35 mm)
Power Requirements	AC 220V/50Hz
Environmental Conditions	0 to 50°C, relative humidity < 80%
Dimensions	230(L)×180(W)×120(H)mm
Weight	2.2 kg

