

Mini TopPette

Pipette

NEW

Features

- **Dual Stop Design:** Enables precise aspiration and dispensing, meeting diverse experimental needs.
- **Accuracy Compliant with ISO 8655:** Ensures reliability and repeatability, enhancing the accuracy of experimental data.
- **Equipped with Tip Ejector:** Allows one-handed operation for easy tip disposal, improving lab efficiency.
- **Color-Coded Plunger Caps by Volume:** Offers intuitive identification of different volume specifications, reducing handling errors.
- **Compact and delicate:** About 13 cm in length, combining comfortable grip with efficient operation.

Specifications

Volume Range	Measured Volume	Systematic Error		Random Error	
	µL	µL	%	µL	%
5	5	±0.125	±2.5	±0.075	±1.5
10	10	±0.12	±1.2	±0.08	±0.8
50	50	±0.5	±1.0	±0.25	±0.5
100	100	±0.8	±0.8	±0.3	±0.3
200	200	±1.6	±0.8	±0.6	±0.3
300	300	±2.4	±0.8	±0.9	±0.3

* For user calibration, please refer to international standard ISO 8655-2.