

显微镜机架 THE STAND OF MICROSCOPES



◀ 20 型机架
20 Series Stand



◀ LBXW 型机架
LBXW Series Stand



◀ LBXY 型机架
LBXY Series Stand



◀ DSM 型机架
DSM Series Stand



▲ DSMW 型机架
DSMW Series Stand



▲ LDX 型机架
LDX Series Stand



▲ TD - I 大平台支架
TD - I Universal Stand With Large Base



▲ TD - II 大平台支架
TD - II Universal Stand With Large Base



▲ TD - III 万向支架
TD - III Universal Stand



▲ TD - IV 万向支架
TD - IV Universal Stand

显微镜配件 SPARE PART FOR MICROSCOPES



▲ 目 镜 Eyepiece




▲ 物 镜 Objective



▲ 附加物镜
Auxiliary objectives




▲ 插入式物镜
Inserted objectives





▲ 载物台
Mechanical stage

显微镜配件
SPARE PART FOR MICROSCOPES






TECHNICAL PARAMETERS 技术参数

型 号 Model	单支光纤 Single Fiber	双支光纤 Dual Fiber	环形光纤 Ring Fiber	
长度 Length(mm)	550	550	900	800
光纤直径 Optical Fiber D(mm)	5.6	5.6	10	8
环形内径 Ring ID(mm)	/	/	55.5	55.5





荧光环形灯
 FLUORESCENT RING LIGHT



型 号 Model	RFL8	RFL8-1	RFL8-2	RFL10-1	RFL10-2	RFL10-2T
电 压 Voltage	110V / 220V					
功 率 Power	8W			10W		
色 温 Color Temperature	2700K / 6000K					
调 光 Light adjustment	无调光 No Light adjustment					调光 Light adjustment

LED 荧光环形灯
LED RING LIGHT



TECHNICAL DATA FOR LED RING LIGHT LED环形灯主要技术参数

[illegible]