











## Optional configuration introduction

Optional configuration	Optional configuration introduction	Picture	Price
Color touch paperless recorder	<ol style="list-style-type: none"> <li>1. 3 – inch LED display, with a standard USB2.0 interface.</li> <li>2. Support USB data transfer function.</li> <li>3. Easy for user to save the data that the GMP certification needed.</li> </ol>		5500RMB
Wireless text messages, alarm system	<ol style="list-style-type: none"> <li>1. If the device goes wrong and user are not in the scene, the system will gather failure promptly to send the text message to the receiving people on the mobile phone, ensure to troubleshoot, recover experiment, avoid damage.</li> <li>2. Can immediately check the working status of equipment and control.</li> <li>3. Excluding mobile phone.</li> </ol>		2000RMB
Embedded, thermal printer	Print the required data		1500RMB
Embedded, needle type printer	<ol style="list-style-type: none"> <li>1. Print the required data</li> <li>2. Printer + controller + power</li> </ol>		1800RMB
16 channel temperature inspecting meter introduction	<ol style="list-style-type: none"> <li>1. Used for temperature uniformity test in the box-type chamber.</li> <li>2. High-precision, 3 meter long, PT100 sensor. Measuring range: 0-300°C</li> <li>3. Backlight type 3.5 inch 128*64 high resolution LCD display screen, 16 channel data display on one screen, can search historical data, change time scale on curve.</li> <li>4. With multiple protection and isolation design, the anti-interference ability is strong and the reliability is high</li> <li>5. Measurement accuracy: 0.2%FS ± 1d, refresh cycle 2S/16 channel.</li> </ol>		4500RMB
Multi – segment programmable LCD temperature controller	<ol style="list-style-type: none"> <li>1. Micro – computer controller, time and heating rate.</li> <li>2. Can preset 30 – step programmable in paragraph 15, each setting time 1-99 hours, 99 minutes.</li> <li>3. Multi – segment programmable control, you can simplify the complex process of experiment, the actual realization of automatic control and operation.</li> </ol>		1500RMB
Liquid rose – independent temperature limit controller	Exceeding the limit temperature interrupt automatically, ensures experiments to run safely.		180RMB
RS485/232 Communication Interface	Display, monitor the data of the experiment via connecting computers.		600RMB
U-Disk data storage, with USB module	USB data transfer interfaces, available with data transfer, with USB module.		1000RMB
BOD socket	BOD instrument – specific socket		150RMB