

# KD-202A Rotary Microtome



## Features

- Easy-to-clean housing
- Blade holder can be moved forward-or-backward and left-or-right resulting in easy trimming and sectioning operations
- LCD screen shows the number of sections and trims; equipped with a safety alert system

## Technical Specifications

Section thickness range	0 – 25 $\mu$ m for continuous sectioning
The sectioning thickness can be set at any value	> 25 $\mu$ m
Minimum Setting of Sectioning Thickness	1 $\mu$ m
Maximum Specimen Size	40×30mm
Precision Error	±5%
Dimensions	340×335×270mm (W×D×H)
Net weight	18.5kg