

# ***JX-5000 CENTRIFUGE***



## **1. GENERAL**

The frame of this instrument adopts electrostatic sprayed plastic. Its model is beautiful, and it has the advantages of small cubage, low weight, large capacity, low noise, high efficiency and so on. It can be used in hospitals and biochemical labs for qualitative analysis of serum, urea and plasma.

## **2. Technical data**

Model	JX-5000
Motor type	DC motor
Max. speed /RCF	Program 1: 5000rpm/3074Xg Program 2: 4000rpm/2220Xg
Display	LCD for all parameters
Capacity	12X15ml, or 8X20ml, or 6X20ml,
Motor power	90W
Power supply	220V 50Hz
Noise	<65dB
Time range	1s-99m59s
Overall size	400X295X225(mm)