

# Shaker

## >>> KJ-201ZT Multifunction Shaker

### Features:

- Reasonable design, the machine will be stable, No shifting or noise while running.
- Adjustable Analog dial design can ensure stepless speed adjustment.
- High performance, low noise DC motor to ensure rapid and thorough mixing.
- Customized plate available to meet different test requirements.



KJ-201ZT

### Technical Data:

Operation method:	Rapid orbital motion
Eccentric diameter:	4mm
Rating speed:	200-2800rpm
Running mode:	Continuous
Max. loading:	15 pcs 50mL centrifuge tubes or 50 pcs 3-10mL tubes
Input voltage(frequency):	AC 220V±10%(50/60Hz)
Permissible relative moisture:	80%
Output power:	10W
Protective class:	IP21
Permissible ambient temperature:	5-50°C
Dimensions:	260×200×260mm
Net. Weight:	2.2Kg