

# 电热套

## Heating Mantle



智能控温  
Intelligent temperature control



加热均匀，升温快  
Heating evenly, heating up fast

### 技术参数

#### TECHNICAL DATA

型号 Model	HM		HMC		HMS			HMSC	
容量 Capacity	100ml	250ml	500ml	1L	2L	3L	5L	10L	20L
功能 Function	加热 Heating		加热 Heating		加热/磁力搅拌 Heating/magnetic stirrer			加热/磁力搅拌 Heating/magnetic stirrer	
加热功率 Heating power	100W	150W	250W	350W	450W	600W	800W	1200W	2400W
控温方式 Temperature control	/		数显控温 Digital temperature control			/		数显控温 Digital temperature control	
加热温度 Max temperature	室温 - 380°C Room temp - 380°C		室温 - 300°C Room temp - 300°C		室温 - 380°C Room temp - 380°C			室温 - 300°C Room temp - 300°C	