

迷你金属浴 Mini Dry Bath



体积小
Small volume



多点运行
Multi-point operation



LCD显示
LCD display



超温保护
Overtemperature protection

技术参数

TECHNICAL DATA

型号 Model	MDB100C	MDB100
温度范围 Temperature range	室温-15°C ~ 100°C Room temp - 15°C ~ 100°C	室温+5°C ~ 100°C Room temp + 5°C ~ 100°C
温度稳定性 Temperature stability	±0.5°C	±0.5°C
温度均匀性 Temperature homogeneity	±0.5°C	±0.5°C
显示精度 Display accuracy	0.1°C	0.1°C
加热时间 Heating time	≤12min (25-100°C)	≤12min (25-100°C)
冷却时间 Cooling time	20 - 25 min (20-0°C)	/
多点运行 Multi-point operation	是 Yes	No
可设温度点 Temperature point can be set	9	/
定时范围 Time setting range	1s - 23h59min	1s - 23h59min
功率 Power	60W	48W
电压 Voltage	AC100-240V	AC100-240V
尺寸 Size	160*115*96 mm	160*115*96 mm

可选模块 Module options	试管数及孔径 Numbers and aperture of test tube
MDB1	40 * 0.2ml(离心管 / Centrifuge tube)
MDB2	24 * 0.5ml(离心管 / Centrifuge tube)
MDB3	15 * 1.5ml(离心管 / Centrifuge tube)
MDB4	15 * 2.0ml(离心管 / Centrifuge tube)
定制模块 OEM of modules	模块可根据客户需求定制 Modules can be customized according to customer needs