

P 系列整只消毒移液器

P Series Whole Autoclavable Micropipette



121°C消毒
121°C



量程范围广
Wide range



分液精准
Dispense accurately



技术参数

TECHNICAL DATA

型号 Model	容量范围 Capacity Range	测试容量 Test Capacity	吸头 Tip	增量 Increment	精度误差 Accuracy Error	重复误差 Repeat Error
P10	0.5 - 10ul	5ul	10ul	0.1ul	±1.00%	±1.00%
P20	2 - 20ul	10ul	200ul	0.5ul	±1.00%	±1.00%
P50	5 - 50ul	25ul	200ul	0.5ul	±0.90%	±0.60%
P100	10 - 100ul	50ul	200ul	1ul	±1.00%	±0.40%
P200	20 - 200ul	100ul	200ul	1ul	±0.80%	±0.30%
P1000	100 - 1000ul	500ul	1000ul	5ul	±0.70%	±0.25%
P5000	1000 - 5000ul	2500ul	5000ul	50ul	±0.60%	±0.30%
P10000	2 - 10ml	5ml	10ml	0.1ml	±1.20%	±0.30%