

BJ-1 DISINTEGRATION TESTER

●BJ series of instruments are used for detecting disintegration of solid in prescriptive conditions.

Applicable standards:

- National standard (Chinese codex) /2015
- Industry Standard (Disintegration tester)

Features:

- There are six tubes in one basket assembly.
- Auto stop at timing point.
- The bath liquid can achieve an even temperature.
- A MPU is used to automatically control the temperature, time and frequency of nacelle back and forth.
- The time parameters can be preset at anytime. The preset and real time data can be displayed alternately in time.
- Automation: auto-test, auto-diagnose, auto-alarm.

Specifications:

- | | |
|--------------------------------------|--------------------------------------|
| ●Basket assembly | 1 |
| ●Nacelle back and forth frequency | (30~32) / min |
| ●Nacelle back and forth range | (55±1) mm |
| ●Automatic time range | (1~900) min |
| ●Temperature range | (ambient to 45) °C |
| ●Accuracy of controlling temperature | (37±0.5) °C |
| ●Power | 220V/ 50Hz/ 600W or 110V/ 60Hz/ 600W |
| ●Dimension | (260×300×420)mm ³ |

