

CM SERIES CLARIFY TESTER

- They test clarify of injection and bottled medicine liquid.

Applicable standards

- National standard(Chinese codex)/2020
- Health department standard
- Corporate standard(Clarify tester)


CM-1 CLARIFY TESTER

Features

- The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor and shading.
- Flash eliminated. Illumination is adjustable.
- Time and alarm automatically. Test continuously.
- Used by one person.

Specifications

- Lamp three-primary-color fluorescent lamp
- Illumination range (1000~4000) LUX
- Power 220V/ 50Hz/ 30W
- Time range (1~99) s



- Dimension (66*18*50) cm

CM-2 CLARIFY TESTER



Features

- The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor and shading.
- Flash eliminated. Illumination is adjustable.
- Time and alarm automatically. Test continuously.
- Used by two person at the same time.

Specifications

- Lamp three-primary-color fluorescent lamp
- Illumination range (1000~4000) LUX
- Time range (1~99) s
- Power 220V/ 50Hz/ 50W
- Dimension (66*38*50) cm