

# **H D — 1 A THICKNESS TESTER**

- HD series of instruments are used for detecting thickness of tablet and capsule.

## **Applicable standard**

- Corporate standard (Thickness tester)

## **Features**

- Used for measuring thickness of peak and wall of capsule.
- High-accuracy test with the most precise, high-quality electronic indicator.

## **Specifications**

- Test range            0~0.5"/12.7 mm
- Test accuracy        0.001"/0.02 mm
- Type of capsule      00、0、1、2、3、4、5 type capsule