

H D — 2 A THICKNESS TESTER

- HD series of instruments are used for detecting thickness of tablet and capsule.

Applicable standard

- Corporate standard (Thickness tester)

Features

- Used for measuring thickness of peak and wall of capsule & diameter of tablet and pill.
- High-accuracy test with the most precise, high-quality electronic indicator.

Specifications

- Test range 0~1"/25.4mm
- Test accuracy 0.001"/0.02 mm
- Type of capsule 00、0、1、2、3、4、5 type capsule
- Dimension (200 × 80 × 50) mm
- weight 0.9 kg

