

TM-1 TRANSPARENCY TESTER

- Transparency tester is requisite for detecting transparency of gelation.

Applicable standards:

- National standard (Food additive Gelatin)
- National standard (Hard gelatin capsule for medicine)
- Vocation standard (Gelatin for medicine)
- Corporate standard (Thickness tester)

Features:

- Thoroughly eliminating the influence to results by environment and experimenter's visual error.
- Standard test tube and same capacity of sample.
- High stable light source system, equal brightness.
- High precise photoelectric transducers for testing exactly.
- Auto-testing, transparency being conversion accounted automatically.
- Automation: auto-test, auto-diagnose, auto-alarm.

Specifications:

- Test range (20~500) mm
- Test accuracy ± 5 mm
- Volume 400ml
- Power 220V/ 50Hz/ 80W ● Dimension (360×240×120) mm³

