

YD-1A TABLET HARDNESS TESTER

- Tablet hardness testers are instruments for detecting breaking hardness of tablet.

Applicable standards:

- Corporate standard (Tablet hardness tester)

Features:

- The tablet is loaded and pressed artificially.
- There is high accuracy press sensor and AD converter.
Hardness data is displayed on LED.
- Automatically zero after warming up.
- Digital calibration and zero adjustment.
- N / KG conversion by one-button.
- It tests tablet hardness continuously.
- Test data is displayed and locked automatically, so can test be reset and cycled.
- Automation: auto-test, auto-diagnose, auto-alarm.

Specifications:

- Hardness range (2~200) N/(0.2~20) kg
- Hardness accuracy $\pm 1.5\%$
- Probe procession range (3~20) mm

- Power 220V/ 50Hz/ 5W or 110V/ 60Hz/ 5W
- Dimension (32*18*11) cm³
- Weight 3.6 kg

