

SEED TESTING INSTRUMENT

High Accuracy Automatic Seed Counter


Model: SLY-D / SLY-E

●Characteristic

- ▲ Microcomputer automatic control system. LCD display; light touch key.
- ▲ It is suitable for measuring round and long, small and big seeds
- ▲ Direct read data on screen, such as current seeds counter number, counting mode, sensitivity, weight and counting time etc.
- ▲ Adjustable counting speed
- ▲ Sensitivity adjustable, make the results more accuracy
- ▲ Two modes for choosing: unlimited counting, and counting stop.
- ▲ Auto save your previous counting
- ▲ Shortcut menu operation "500", "1000", "1500", "2000" options.
- ▲ 4G SD card including.
- ▲ Direct read data on screen, such as count grain number, date, time, time consuming, weight
- ▲ The weighing function. SLY-E type can automatic weighing.
- ▲ Can display and save the results

●Technical specification

Counting range: L*W less than 12*4mm (small seeds)
L*W less than 12*10mm (big seeds)
Counting accuracy: 5/1000 for small seeds; 3/1000 for big seeds
Counting speed: 1000pcs/ 3m
Counting capacity: 1-99999 (read through LCD)
Preset self stop: any value within 1-99999
Weighing range: 0-1000g
Size: 450*300*320mm
Power supply: AC220±20V~ 50Hz
Watt: <60W
Continuous working time: ≥ 5h



● Model Difference

| Model | Function |
|-------|---|
| SLY-D | with SD card, data can be stored and export. without the weighing function. |
| SLY-E | with SD card, data can be stored and export. with the weighing function. |